

PULSED AF¹ Pivotal Trial Overview

The PulseSelect™ pulsed field ablation (PFA) system is limited to investigational use only and is not approved for commercial use in any geography. CE Mark pending.

Trial design and study population

Trial design

Paired single arm. Prospective, nonrandomised clinical study.

Global multicentre study

9 countries: Austria, Belgium, France, Spain, Netherlands, United States, Canada, Australia, Japan.

41 sites

67 operators

300 subjects



150 Paroxysmal (PAF)

150 Persistent (PsAF)


300
patients

The study population included recurrent symptomatic paroxysmal and persistent atrial fibrillation patients refractory to class I or III antiarrhythmic drugs.

96%

96% (287) of patients completed 12-month follow-up.

Rigorous arrhythmia monitoring



Acute procedural failure modes

Failure modes

Trial outcomes

PULSED AF evaluated the safety and effectiveness of the PulseSelect™ Pulsed Field Ablation System for the treatment of patients with paroxysmal or persistent atrial fibrillation.



0 Esophageal events



0 PV stenosis



0 Phrenic nerve injury



0 Coronary artery spasm

**Primary
safety endpoint**

0.7%

Primary adverse
events (95% CI, 0.1
to 4.6%)

Performance goal of < 13%

1/300 Cerebrovascular accident

1/300 Tamponade

0/300 Transient ischemic attack

0/300 Major bleeding

0/300 Myocardial infarction

0/300 Pericarditis

0/300 Vagal nerve injury

0/300 Systemic pulmonary embolism

0/300 Pulmonary edema

0/300 Vascular access complications

0/300 Cardiovascular hospitalisation

0/300 Death

Quality-of-life scores improved post-ablation compared to baseline




AFEQT score improved by 29.4 (95% CI, 25.8 to 33.1) and 29.0 (95% CI, 25.5 to 32.5) points in the paroxysmal and persistent populations respectively from baseline to 12 months.



EQ-5D-5L score improved by 0.05 (95% CI, 0.02 to 0.08) points in paroxysmal **and 0.06 (95% CI, 0.04 to 0.09) points** in persistent atrial fibrillation patients.



12-month effectiveness data

Freedom from any atrial arrhythmias recurrence of ≥ 30 seconds



Primary effectiveness

Acute procedure failure, AF/AFL/AT recurrence, cardioversion, repeat ablation, new/re-initiated/increased AADs, any subsequent AF surgery




Clinical success


Freedom from recurrence of any symptomatic atrial arrhythmias (post-hoc analysis)



Procedural data (ancillary endpoints)



Procedure times under 50 minutes
when excluding the 20-minute trial-mandated wait period.



Total PFA energy delivery under 30 seconds

Parameter	Paroxysmal (n = 150)	Persistent (n = 150)
Skin-to-skin procedural time (min) [†]	134 ± 50	145 ± 60
Device left atrial dwell time (min) [‡]	65 ± 29	70 ± 31
Fluoroscopy time during procedure (min)	26 ± 17 [§]	29 ± 21

The PulseSelect™ system was used in conjunction with multiple commercially available mapping systems.

Plus-minus values are means ±SD, with median (interquartile range) also presented.


[†]Skin-to-skin procedure time is from first sheath inserted to last sheath pulled out – this includes transfer time to recovery room and sheath pulling time.

[‡]Left atrial dwell time and time from first to last application includes the protocol-mandated 20-minute wait period and any post-ablation mapping time.

[§]Data were available for 149 patients.

Stepwise ostial-antral approach

- Four pulse trains delivered per application
- One application per catheter position
- A minimum of 8 applications per vein is recommended to achieve pulmonary vein isolation.



¹ Verma A, Haines DE, Boersma LV, et al. Pulsed Field Ablation for the Treatment of Atrial Fibrillation: PULSED AF Pivotal Trial. *Circulation*. May 9, 2023;147(19):1422-1432.

Medtronic

Medtronic International Trading Sarl
Route du Molliau 31
Case postale
1131 Tolochenaz
Switzerland
Tel: +41 (0) 21 802 70 00
Fax: +41 (0) 21 802 79 00

United Kingdom/Ireland
Medtronic Limited
Building 9
Croxley Park
Hatters Lane
Watford
Herts WD18 8WW
www.medtronic.co.uk
Tel: +44 (0)1923 212213
Fax: +44 (0)1923 241004

medtronic.eu

The products are limited to investigational use only and are not approved for commercial use in any geography. CE Mark pending. ©2023 Medtronic. Medtronic, Medtronic logo, and Engineering the extraordinary are trademarks of Medtronic. All other brands are trademarks of a Medtronic company. UC202404559-pulsed-af-pivotal-trial-overview-brochure-en-emea-11019703