

INTELLANAV™ ST

ABLATION CATHETER

Familiarity
Proven Performance
Confident Navigation



Precision of a small solid tip with magnetic navigation accuracy for precise ablations in SVT procedures, critical when near the conduction system.

RHYTHMIA HDx™ MAPPING SYSTEM Accessories

Model Number	Description
M004 RARC01 0	INTELLANAV Ablation Catheter Cable 2m
M004 RA6201 0	RHYTHMIA HDx Connection Box - MAESTRO

RHYTHMIA™ MAPPING SYSTEM Accessories

Model Number	Description
M004 RARC01 0	INTELLANAV Ablation Catheter Cable 2m
M004 RA6200 0	INTELLANAV Connection Box - MAESTRO
M004 RARC20 0	INTELLANAV SIU Adapter Cable
M004 5441S 0	Octapolar Dx Cable



Ordering Information: **INTELLANAV™ ST** Electrode Configuration: Quadripolar Electrode Spacing: 2.5 mm

Model Number	Shaft Size	Tip Size	Curve Style	Shaft Length
M004 EPTR5031TH 0*	7F	7F/4mm	Standard	110cm
M004 EPTR5031THK2 0*	7F	7F/4mm	Large	110cm

7F = 2.33mm

Bidirectional Curve Options



Boston Scientific curve Substitutions

BSC Standard	Biosense: D Curve Abbott: D Curve, Medium Curl, Medium Sweep
BSC Large	Biosense: F Curve, J Curve Abbott: F Curve, J Curve, Large Curl, Large Sweep, Large

Catheter configurations are illustrative representations only and may not reflect actual performance.
*Internal diameter dimensions: Standard: 25mm / Large: 33mm

All cited trademarks are the property of their respective owners. CAUTION: The law restricts these devices to sale by or on the order of a physician. Indications, contraindications, warnings and instructions for use can be found in the product labeling supplied with each device. Information for the use only in countries with applicable health authority product registrations. Material not for distribution and use in France. EP-572609-AA © 2018 Boston Scientific Corporation or its affiliates. All rights reserved. Printed in Germany by medicalvision. DINEP2292EA